



Electrothermal Electric Bunsen Heating Mantles

Electrothermal Electric Bunsen heating mantles are corrosion-resistant and include easy-to-replace heating elements.

Combines the advantages of a gas burner with the clean operation and easy control of electric heating.

Radiation from the heater is directed upwards only, so the Bunsen burner is ideal for heating test tubes, crucibles, small flasks and beakers, independent of their shape.

- Burner consumes only 400w of power
- Conical shaped heating element directs radiant heat to the top cavity
- Top cowl deflects heat away from your hand
- Air circulation from the vented housing keeps the base cool enough to hold during operation

Warranty: One year parts and labour

Certifications: CE

Specifications	
Capacity	25ml (0.006gal)
Material	Stainless-Steel
Temperature Range	800° to 1000° C
Flask Volume	25ml
Shipping Weight	.975kg (2.2lbs)
Dia x H	12 x17.7cm (4.5 x7in.)

Without Controller	
BA6101	230V 50/60Hz, 480w
BA6101X1	115v 50/60Hz, 430w
BA6101X2	100v 50/60Hz, 480w
BA6101X3	90v 50/60Hz, 430w
BA6101X6	230v 50/60Hz, 480w, EU Plug
With Controller	
BA6101/C	230v 50/60Hz, 480w
BA6101/CX1	115v 50/60Hz, 430w
BA6101/CX6	230v 50/60Hz, 480w, EU Plug