



Electrothermal Paraffin Section Mounting Bath

The Electrothermal Paraffin Section Mounting Bath allows for comfortable and easy manipulation of sections.

- The PTFE black interior and stove enamelled white aluminium exterior is virtually maintenance free and extremely easy to clean.
- Excess or damaged sections may be simply removed by skimming a filter paper over the water surface

Operational temperature and performance.

- 70 °C Nominal with 2 litres water. Obtained within 120 minutes.

Warranty: Once year parts and labour

Certifications: CE

Specifications		
Material	Aluminium base with PTFE Xylan coated interior.	
Temperature	Ambient +5 °C to 75 °C using arbitrary 1-10 scale.	
Cat. No.	Capacity	Electrical Requirement
MH8516	2.25 litres (0.6 gallon)	230V 50/60Hz
MH8516/X1	2.25 litres (0.6 gallon)	115v 50/60Hz
MH8516/X6	2.25 litres (0.6 gallon)	230v 50/60Hz, EU Plug

Lid Sold Separate (Cat. No. AZ9241)