

Overhead Stirrers

Stands

Extra strong retort stands with H-pattern base for stability and robust support rod. Both versions include bosshead. Heavy duty stand recommended for use with SS30 overhead stirrer.

Technical Specification

	SS10/1	SS10/2
Base, mm (w x d x h)	400 x 350 x 25	550 x 480 x 25
Rod, mm (dia. x l)	16 x 700	25 x 850
Net weight, kg	7.6	20.2

Ordering Information

Model	Description
SS10/1	Stand, general purpose
SS10/2	Stand, heavy duty

Stirrer paddles

Made from high grade stainless steel, there is a choice of five paddle heads and two lengths of 8mm diameter paddle rod. The heads screw on and off the rods, so they can be mixed and matched depending upon stirring requirements. This gives maximum versatility and value for money.

For complete paddle, order rod plus head(s).

Technical Specification

Paddle rods

	SS10/5	SS10/6
Diameter	8mm	8mm
Length	350mm	550mm

Paddle heads

	SS10/10	SS10/11	SS10/12	SS10/13	SS10/14
Total width	60mm	94mm	80mm	50mm	60mm

Ordering Information

Model	Description
SS10/5	Paddle rod, 350mm
SS10/6	Paddle rod, 550mm
SS10/10	Small paddle head
SS10/11	Large paddle head
SS10/12	Large propeller head (2 blade)
SS10/13	Small propeller head (4 blade)
SS10/14	Anchor paddle head

