

Electrothermal Slow Speed Stirrers with Integral Control



Electrothermal Slow Speed Stirrers are ideal for use with magnetic culture vessels.

Stirrers feature sophisticated, electronically controlled magnetic drivers for precision stirring, which is reproducible to better than 1 rpm.

- Splash-proof noncorrosive polymer case
- Low-profile, easy-to-clean design

Warranty: Three years

Certifications: CE

Specifications	
Stirring Range	0 to 150 RPM

Cat. No.	No. of Stirring Positions	Capacity	Overall L x W x H	Weight	Electrical Requirement/ Plug Type
PS60046	1	2L	15x15x7cm (6x6x2.75in)	0.8kg (1.8lb.)	230V/UK
PS60046X1	1	2L	15x15x7cm (6x6x2.75in)	0.8kg (1.8lb.)	115V/US
PS60046X6	1	2L	15x15x7cm (6x6x2.75in)	0.8kg (1.8lb)	230V/EU
PS60058	4	1L per position	26x26x7cm (10x10x2.75in)	2.7kg (6lb)	230V/UK
PS60058X1	4	1L per position	26x26x7cm (10x10x2.75in)	2.7kg (6lb)	115V/US
PS60058X6	4	1L per position	26x26x7cm (10x10x2.75in)	2.7kg (6lb)	230V/EU
PS60060	1	10L	26x26x7cm (10x10x2.75in)	2kg (4.4lb)	230V/UK
PS60060X1	1	10L	26x26x7cm (10x10x2.75in)	2kg (4.4lb)	115V/US
PS60060X6	1	10L	26x26x7cm (10x10x2.75in)	2kg (4.4lb)	230V/EU
PS60055	1	24L	26x26x7cm (10x10x2.75in)	2kg (4.4lb)	230V/UK
PS60055X1	1	24L	26x26x7cm (10x10x2.75in)	2kg (4.4lb)	115V/US
PS60055X6	1	24L	26x26x7cm (10x10x2.75in)	2kg (4.4lb)	230V/EU