



SA7

Mixer, vortex, fixed speed, SA7

- Fixed speed 2500rpm
- Robust die-cast construction
- Automatic press start

For rapid mixing of samples contained in test-tubes, small flasks and bottles. Starts automatically when the rubber cup is depressed and stops once vessel is removed.

Heavy die-cast body for stability during use.

Technical Specification

	SA7
Speed	Fixed 2500rpm
Orbit diameter, mm	4.2
Dimensions, mm (w x d x h)	135 x 215 x 78
Net weight, kg	3.2
Electrical supply	90 – 240V, 50 / 60Hz, 20W
IP Rating	31

Ordering Information

Model	Description
SA7	Vortex mixer, fixed speed